



Configuration du MACS Quant 10



Configuration du MACS Quant 10

489→ **GFP** ↗510
500→ **Lysotracker Green** ↗510
493→ **Calcein** ↗514
495→ **FITC** ↗519
499→ **Alexa Fluor 488** ↗520
495→ **CFSE** ↗525
496→ **Vio Bright FITC** ↗522
490→ **BD Horizon BB515** ↗515
490→ **Zombie Green** ↗515
498→ **BD FVS 520** ↗521
514→ **YFP** ↗527

