


Configuration du LSR FORTESSA (PT Cytométrie Purpan)

Blue Laser (488 nm)

↻ 477 **PerCP** 678 ↻
↻ 477 **PerCP-Cy5.5** (↻ 675) 696 ↻
↻ 477 **PerCP-Vio700** 704 ↻
↻ 477 **PerCP-eFluo710** 710 ↻

SSC



↻ 489 **GFP** 510 ↻
↻ 500 **Lysotracker Green** 510 ↻
↻ 493 **Calcein** 514 ↻
↻ 495 **FITC** 519 ↻
↻ 499 **Alexa Fluor 488** 520 ↻
↻ 495 **CFSE** 525 ↻
↻ 496 **VioBright FITC** 522 ↻
↻ 490 **BD Horizon BB515** 515 ↻
↻ 490 **Zombie Green** 515 ↻
↻ 498 **BD FVS 520** 521 ↻

Configuration du LSR FORTESSA


(PT Cytométrie Purpan)

Yellow Green Laser

(561 nm)

↻ 565 **PE-Vio615** 519 ↻
 ↻ 565 **PE-eFluor610** 607 ↻
 ↻ 566 **PE-Dazzle594** 610 ↻
 ↻ 496 **PE-CF594** 612 ↻
 ↻ 566 **PE-Texas-Red (ECD)**(↻ 595) 616 ↻
 ↻ 566 **PE-Alexa Fluor 610** (↻ 610) 627 ↻
 ↻ 351 **PI** 617 ↻

↻ 566 **PE-Cy5.5** (↻ 675) 693 ↻



↻ 547 **Cell Tracker Orange (CMTMR)** 575 ↻

↻ 566 **PE** 576 ↻

↻ 553 **Ds-Red** 586 ↻

↻ 565 **PE-Vio770** 775 ↻

↻ 566 **PE-Cy7** 778 ↻

↻ 543 **7-AAD** 647 ↻

↻ 566 **PE-Cy5** (↻ 647) 666 ↻


Configuration du LSR FORTESSA (PT Cytométrie Purpan)

Violet Laser (405 nm)

↪ 300 **Qdot 705** 702 ↪
↪ 405 **BV711** 711 ↪

↪ 300 **Qdot 605** 603 ↪
↪ 405 **BV605** 603 ↪
↪ 405 **eVolve605** 610 ↪

↪ 405 **BV421** 421 ↪
↪ 400 **Zombie Violet FV423** ↪
↪ 398 **Cascade Blue** 420 ↪
↪ 401 **Alexa Fluor 405** 422 ↪
↪ 404 **Pacific Blue** 456 ↪
↪ 360 **Dapi Viab** 461 ↪
↪ 400 **VioBlue** 452 ↪
↪ 405 **eFluor450** 450 ↪
↪ 404 **BD Horizon V450** 448 ↪
↪ 406 **BD FVS 450** 450 ↪



↪ 458 **AmCyan** 489 ↪
↪ 400 **Cascade Yellow** 550 ↪
↪ 400 **Pacific Orange** 551 ↪
↪ 405 **BV510** 510 ↪
↪ 415 **BD Horizon V500** 500 ↪
↪ 388 **VioGreen** 520 ↪
↪ 408 **BD FVS 510** 512 ↪
↪ 385 **Zombie Aqua FV** 515 ↪

↪ 300 **Qdot 655** 654 ↪
↪ 405 **BV650** 645 ↪
↪ 405 **eVolve655** 655 ↪


↪ 300 **Qdot 800** 792 ↪
↪ 405 **BV785** 785 ↪

Configuration du LSR FORTESSA (PT Cytométrie Purpan)

Red Laser (640 nm)

↪ 650 **APC-Alexa Fluor 750** (↪ 752) 779↪
↪ 650 **APC-H7** 780↪
↪ 650 **APC-Cy7** (↪ 753) 785↪
↪ 749 **Alexa Fluor 750** 775↪
↪ 652 **APC-Vio770** 775↪
↪ 650 **APC-eFluor780** 780↪
↪ 759 **BD FVS780** 780↪

↪ 650 **APC** 660↪
↪ 631 **Alexa Fluor 633** 647↪
↪ 630 **eFluor660** 660↪
↪ 650 **Alexa Fluor 647** 668↪
↪ 647 **Cy5** 660↪
↪ 649 **BD FVS 660** 660↪



↪ 650 **APC-Cy5.5** (↪ 675) 695↪
↪ 675 **Cy5.5** 695↪
↪ 680 **Alexa Fluor 680** 700↪
↪ 700 **Alexa Fluor 700** 720↪
↪ 652 **APC-R700** 704↪
↪ 657 **BD FVS700** 700↪