


Configuration du FUSION

Violet Laser (405nm)

407 BV711 711
414 Super Bright 702 702

407 BV605 605
414 Super Bright 600 600
405 eVolve 605 610

407 BV421 421
400 DyLight 405 420
402 Alexa Fluor 405 421
414 Super Bright 436 436
404 BD Horizon V450 448
405 eFluor 450 450
400 VioBlue 452
404 Pacific Blue 456
CTV



436 BV480 478
405 BV510 510
415 BD Horizon V500 500
388 VioGreen 520
431 Alexa Fluor 430 541
400 Pacific Orange 551

407 BV650 650
414 Super Bright 645 645
405 eVolve 655 655


407 BV786 786
407 BV785 786
414 Super Bright 780 780

Configuration du FUSION

Red laser (637nm)

650 APC-H7 780
650 APC-Cy7 (753) 785
650 APC-Alexa Fluor 750 (752) 779
749 Alexa Fluor 750 775
652 APC-Vio770 775
650 APC-eFluor780 780
759 BD FVS780 780

650 APC 660
631 Alexa Fluor 633 647
630 eFluor660 660
633 Alexa Fluor 635 647
650 Alexa Fluor 647 668
647 Cy5 660
649 BD FVS 660 660
654 DyLight 650 673




652 APC-R700 704
679 Alexa Fluor 680 702
692 DyLight 680 712
702 Alexa Fluor 700 723
650 APC-Cy5.5 695

Configuration du FUSION

Blue Laser (488 nm)

BB700
477 PerCP 678
477 PerCP-Cy5.5 (↗ 675) 696
477 PerCP-Vio700 704
477 PerCP-eFluo710 710

SSC




BB515
489 GFP 510
500 LysoTracker Green 510
493 Calcein 514
495 FITC 519
499 Alexa Fluor 488 520
495 CFSE 525
496 VioBright FITC 522
490 BD Horizon BB515 515
490 Zombie Green 515
498 BD FVS 520 521

Configuration du FUSION

Yellow Green Laser (561nm)

BYG790-P
565 PE-Vio770 775
566 PE-Cy7 778


BYG584-P
566 PE 576
553 Ds-Red 586
547 Cell Tracker Orange 575



566 PE-CF594 612
565 PE-Vio615 519
578 Alexa Fluor 568 603
PE-eFluor® 610 607
566 PE Dazzle 594 610
566 PE-Texas-Red 616
612 Alexa Fluor 610 628
585 mCherry 610

Configuration du FUSION

UV Laser (355nm)



BUV496
382 Zombie Aqua 515

