


Configuration du FUSION

Violet Laser (405nm)




Configuration du FUSION

Red laser (637nm)

650 APC-H7 780
650 APC-Cy7 (753) 785
650 APC-Alexa Fluor 750 (752) 779
749 Alexa Fluor 750 775
652 APC-Vio770 775
650 APC-eFluor780 780
759 BD FVS780 780

650 APC 660
631 Alexa Fluor 633 647
630 eFluor660 660
633 Alexa Fluor 635 647
650 Alexa Fluor 647 668
647 Cy5 660
649 BD FVS 660 660
654 DyLight 650 673




652 APC-R700 704
679 Alexa Fluor 680 702
692 DyLight 680 712
702 Alexa Fluor 700 723
650 APC-Cy5.5 695

Configuration du FUSION

Blue Laser
(488 nm)

BB700
477 PerCP 678
477 PerCP-Cy5.5 (↗ 675) 696
477 PerCP-Vio700 704
477 PerCP-eFluo710 710


SSC



BB515
489 GFP 510
500 LysoTracker Green 510
493 Calcein 514
495 FITC 519
499 Alexa Fluor 488 520
495 CFSE 525
496 VioBright FITC 522
488 VioBright 515 515
488 Vio 515 515
490 BD Horizon BB515 515
490 Zombie Green 515
498 BD FVS 520 521


Configuration du FUSION

Yellow Green Laser (561nm)



Configuration du FUSION

UV Laser
(355nm)



BUV496
382 Zombie Aqua 515

